

CONTOURED

MASSAGE COUCH
2 TO 5 SECTIONS



Ergonomic contoured surface · electric height adjustment · available in 2 to 5 section configurations

The Contoured Massage Couch is built for therapists who demand both clinical performance and patient comfort. Available in multiple configurations from 2 to 5 sections. Electric hi-lo height adjustment (48–93.5 cm) via footswitch. Headrest adjusts -75° to $+30^{\circ}$ via gas strut. Padded face rest and prone arm support included. 70 cm wide platform. 250 kg capacity. TiMotion motor. 7-year frame warranty. CE, FDA, ISO 13485:2016.

 0410 148 292



www.physiomed.au



info@physiomed.au

CONTOURED

MASSAGE COUCH
2 TO 5 SECTIONS



Specification	Details
Model	CM
Configurations Available	2, 3, 4, or 5 sections
Total Length	195 cm
Width	70 cm
Height Range	48 cm – 93.5 cm (electric, ~30 sec full travel)
Headrest Adjustment	-75° to +30° (gas-operated strut)
Padding Thickness	5 cm
Foam Density	36 kg/m ³
Weight Capacity	250 kg
Upholstery	Fire, stain, mildew, oil & wear resistant PVC vinyl
Castors	7.5 cm retractable with safety lock
Motor	TiMotion — electric Hi-Lo via footswitch
Face Rest	Padded face rest included
Prone Arm Support	Included
Frame Legs	Rounded for patient safety
Colours	Navy Blue · Black
Warranty	7 yr frame · 3 yr electrics · 12 mth footswitch
Certifications	CE · FDA · ISO 13485:2016

 0410 148 292



www.physiomed.au



info@physiomed.au