

# FISIOTEK 3000 N

Product Code: F3N | Manufacturer: Rimec, Italy

## Knee & Hip CPM — Built for Clinic and Home Rental



The Fisiotek 3000 N is a Continuous Passive Motion (CPM) device from Rimec, designed to mobilise the knee and hip through controlled flexion-extension movement. Its asymmetric frame provides anatomically correct positioning throughout the range of motion, accommodating patient limbs between 72 and 100 cm.

Specifically noted for its solid, durable construction — making it a preferred choice for home rental programs. Control is available via a stationary console and a mobile programming keyboard. A Warm Up mode gently prepares the joint before therapeutic ranges are applied.

### KEY FEATURES

- Knee and hip mobilisation in a single device
- Solid construction recommended for home rental programs
- Asymmetric frame for anatomically correct positioning
- Dual control: stationary console + mobile patient keyboard
- Warm Up mode for pre-session joint preparation
- Fully programmable: amplitude, velocity, force, duration, pause
- High-quality aluminium and stainless steel construction

# FISIOTEK 3000 N

## TECHNICAL SPECIFICATIONS

Specification	Value
Joints Treated	Knee, Hip
Knee Flexion-Extension	0° to 110°
Hip Range	7° to 115°
Limb Accommodation	72–100 cm (optional accessory to 61 cm)
Power Supply	90–250V ~ 50–60 Hz
Safety Classification	Class I B (EN 60601-1)
EMC	Group 1 Class B (EN 60601-1-2)
Device Classification	CE Class IIa (CEE 93/42)
Weight	14 kg
Programmable Functions	Amplitude, velocity, force, duration, pause, warm up

## CLINICAL INDICATIONS

Knee arthroplasty, meniscectomy, meniscal repair, fracture management, osteotomy, arthrolysis, hip prosthetic surgery.

Supplied exclusively in Australia by Physiomed Australia