

# FIVE SECTION

TREATMENT COUCH  
WITH POSTURAL DRAINAGE



## Motorised postural drainage + five-section positioning — the clinical choice for respiratory physio

The ET52 is designed for physiotherapists who need maximum patient positioning versatility. Electric hi-lo height adjustment (50–96 cm) via footswitch. Motorised postural drainage to +25° for gravity-assisted airway clearance. Five independently adjustable sections — backrest to +70°, headrest –35° to +70° via gas struts. Face hole with plug. Height-adjustable padded armrests. 250 kg capacity. 7-year frame warranty. CE, FDA, ISO 13485:2016 certified.

0410 148 292



[www.physiomed.au](http://www.physiomed.au)



[info@physiomed.au](mailto:info@physiomed.au)

# FIVE SECTION

TREATMENT COUCH  
WITH POSTURAL DRAINAGE



Specification	Details
Model	ET52
Total Length	195 cm
Width	66 cm
Headrest Length	43 cm
Mid Section	51 cm
Backrest Length	96 cm
Backrest Adjustment	0° to +70° (gas strut)
Headrest Adjustment	-35° to +70° (gas strut)
Height Range	50 cm – 96 cm (electric, ~30 sec full travel)
Postural Drainage	+25° (motorised, via footswitch)
Weight Capacity	250 kg
Foam	5 cm / 36 kg/m <sup>3</sup>
Upholstery	Fire, stain, mildew, oil & wear resistant PVC vinyl
Armrests	Height-adjustable padded armrests
Castors	7.5 cm retractable with safety lock
Motor	TiMotion — electric Hi-Lo via footswitch
Face Hole	Yes — with plug
Warranty	7 yr frame · 3 yr electrics · 12 mth footswitch
Certifications	CE · FDA · ISO 13485:2016

 0410 148 292



[www.physiomed.au](http://www.physiomed.au)



[info@physiomed.au](mailto:info@physiomed.au)