

# CLINICAL STOOLS

GAS LIFT  
ROUND & SADDLE



## Ergonomic seating for the treating clinician — round and saddle configurations

Physiomed gas lift stools protect the clinician across long treatment sessions. Round stool: 35 cm seat, 51–72 cm height, optional backrest and foot ring. Saddle stool: 42 cm wide,  $-3^{\circ}$  to  $+6^{\circ}$  seat angle, optional backrest. All models: 160 kg capacity, five-star base, pneumatic gas lift.

Available in Navy Blue, Black, and Grey.

 0410 148 292



[www.physiomed.au](http://www.physiomed.au)



[info@physiomed.au](mailto:info@physiomed.au)

# CLINICAL STOOLS

GAS LIFT  
ROUND & SADDLE



Specification	Details
Round Stool Standard — Height	51 cm – 67 cm
Round Stool Standard — Seat Diameter	35 cm
Round Stool Premium (Foot Ring) — Height	52 cm – 72 cm
Round Stool with Backrest — Height	47 cm – 63 cm
Round Stool Backrest	36.5 x 26.5 cm · 25° angle adjustment
Saddle Stool Standard — Seat	42 cm wide
Saddle Stool — Seat Angle	-3° to +6°
Saddle Stool with Backrest — Height	47 cm – 63 cm
Saddle Stool Backrest	36.5 x 26.5 cm · 25° angle adjustment
Weight Capacity (all models)	160 kg
Lift Mechanism	Pneumatic gas lift
Base	Five-star for stability
Foam	4–6 cm depending on model
Colours	Navy Blue · Black · Grey

 0410 148 292



[www.physiomed.au](http://www.physiomed.au)



[info@physiomed.au](mailto:info@physiomed.au)