



## **120 cm wide — purpose-built for Bobath, NDT, and neurological physiotherapy**

At 120 cm wide — nearly twice a standard table — the NT120 gives the neurological physiotherapist room to work alongside, behind, and across the patient using Bobath and NDT techniques. Three face holes with plugs accommodate all patient heights in prone. Electric height adjustment 49–93.5 cm. Backrest 0°–65° via gas strut. 250 kg capacity. 7-year frame warranty.

# NEUROLOGICAL

BOBATH TREATMENT  
COUCH — NT120



Specification	Details
Total Length	190 cm
Width	120 cm — extra wide Bobath design
Height Range	49 cm – 93.5 cm (electric)
Backrest Angle	0° to +65° (gas-operated strut)
Face Holes	3 — with individual plugs
Weight Capacity	250 kg
Foam	5 cm / 36 kg/m <sup>3</sup>
Upholstery	Fire, stain, mildew, oil & wear resistant PVC vinyl
Castors	7.5 cm retractable
Motor	TiMotion — electric Hi-Lo via foot bar
Colours	Navy Blue
Warranty	7 yr frame · 3 yr electrics · 12 mth footswitch
Certifications	CE · FDA · ISO 13485:2016



0410 148 292



[www.physiomed.au](http://www.physiomed.au)



[info@physiomed.au](mailto:info@physiomed.au)