

STEALTH

PHYSIO TREATMENT
COUCH



Premium black-frame aesthetics — the same proven clinical performance as the ET3

The Stealth is the premium black-frame evolution of the ET3. Short 40 cm head section for unrestricted cervical access. Electric height 47.5–93 cm via 360° foot bar. Back section 0°–80°, head section –45° to +40° via dual gas struts. Face hole with plug. Contoured black PVC vinyl upholstery. 250 kg capacity. TiMotion motor. 7-year frame warranty. CE, FDA, ISO 13485:2016.

 0410 148 292



www.physiomed.au



info@physiomed.au

STEALTH

PHYSIO TREATMENT
COUCH



Specification	Details
Total Length	195 cm
Width	66 cm
Head Section	40 cm (short head)
Mid Section	54 cm
Leg Section	93 cm
Height Range	47.5 cm – 93 cm (electric, ~30 sec full travel)
Back Section Angle	0° to +80° (dual gas struts)
Head Section Angle	-45° to +40° (gas struts)
Weight Capacity	250 kg
Foam	5 cm, 36 kg/m ³ density
Upholstery	Heavy-duty black PVC vinyl — contoured design
Castors	7.5 cm retractable with safety lock
Motor	TiMotion — electric Hi-Lo via 360° foot bar
Frame	Sleek black frame with rounded safety legs
Face Hole	Yes — with plug
Colour	Black (Stealth Series)
Warranty	7 yr frame · 3 yr electrics · 12 mth footswitch
Certifications	CE · FDA · ISO 13485:2016

 0410 148 292



www.physiomed.au



info@physiomed.au