

THREE SECTION

SHORT HEAD PHYSIO TREATMENT TABLE



The first choice for physiotherapy — unrestricted cervical access, full electric height adjustment

The 40 cm short head section gives the treating physiotherapist unrestricted access around the patient's head, neck, and cervical spine — a defining feature of a true physio table. Electric height adjustment (47.5–93 cm) via foot switch protects practitioner OH&S. Back section adjusts 0°–80°, head section –45° to +40° via dual gas struts. Face hole with plug, retractable castors, TiMotion motor. 250 kg capacity. 7-year frame warranty.

 0410 148 292

 www.physiomed.au

 info@physiomed.au

THREE SECTION

SHORT HEAD PHYSIO TREATMENT TABLE



Specification	Details
Total Length	195 cm
Width	66 cm
Head Section	40 cm (short head)
Mid Section	54 cm
Leg Section	93 cm
Height Range	47.5 cm – 93 cm (electric)
Back Section Angle	0° to +80° (dual gas struts)
Head Section Angle	-45° to +40° (gas struts)
Weight Capacity	250 kg
Foam	5 cm / 36 kg/m ³
Upholstery	Fire, stain, mildew, oil & wear resistant PVC vinyl
Castors	7.5 cm retractable
Motor	TiMotion — foot switch / 360° foot bar
Face Hole	Yes — with plug
Colours	Navy Blue · Black · Stealth (black frame)
Warranty	7 yr frame · 3 yr electrics · 12 mth footswitch
Certifications	CE · FDA · ISO 13485:2016

 0410 148 292



www.physiomed.au



info@physiomed.au