

# TILT TABLE

TWO SECTION  
MOTORISED TILT



## Progressive standing rehabilitation — motorised tilt to 85° vertical with full safety system

Essential for neurological, spinal cord injury, and cardiovascular rehabilitation. Motorised tilt to +85° via hand switch. Safety straps on integrated strap rail. Individually adjustable foot plates -10° to +20°. Activity tray included. Electric height 59–104 cm. Head section to +75°. Face hole with plug. 225 kg capacity. 7-year frame warranty.

0410 148 292



[www.physiomed.au](http://www.physiomed.au)



[info@physiomed.au](mailto:info@physiomed.au)

# TILT TABLE

TWO SECTION  
MOTORISED TILT



Specification	Details
Total Length	194 cm
Width	64 cm
Height Range	59 cm – 104 cm (electric, ~30 sec)
Tilt Range	Motorised 0° to +85° (vertical)
Head Section Angle	0° to +75° (gas-operated strut)
Foot Plate Angle	-10° to +20° (individually adjustable)
Weight Capacity	225 kg
Controls	Hand switch — electric Hi-Lo and tilt
Safety Straps	Multiple — mounted on integrated strap rail
Activity Tray	Included
Face Hole	Yes — with plug
Foam	5 cm / 36 kg/m <sup>3</sup>
Upholstery	Fire, stain, mildew, oil & wear resistant PVC vinyl
Castors	7.5 cm retractable
Warranty	7 yr frame · 3 yr electrics · 12 mth footswitch
Certifications	CE · FDA · ISO 13485:2016



 0410 148 292



[www.physiomed.au](http://www.physiomed.au)



[info@physiomed.au](mailto:info@physiomed.au)