

TWO SECTION

ELECTRIC
TREATMENT TABLE



Simple, strong, and built for the rigours of daily clinical use

One of Australia's most popular treatment tables. The two-section design delivers everything a busy physiotherapy or general practice needs: electric height adjustment (47–92 cm), a gas strut head section (-20° to $+75^{\circ}$), face hole with plug, and a 71 cm wide platform for superior lateral support. 250 kg capacity, 7-year frame warranty. Available in Navy Blue, Black, and Grey.

0410 148 292



www.physiomed.au



info@physiomed.au

TWO SECTION

ELECTRIC
TREATMENT TABLE



Specification	Details
Total Length	195 cm
Width	71 cm
Head Section	67 cm
Leg Section	123 cm
Height Range	47 cm – 92 cm (electric)
Head Section Angle	-20° to +75° (dual gas struts)
Weight Capacity	250 kg
Foam	5 cm / 36 kg/m ³
Upholstery	Fire, stain, mildew, oil & wear resistant PVC vinyl
Castors	7.5 cm retractable
Motor	TiMotion — electric Hi-Lo via foot switch
Face Hole	Yes — with plug
Colours	Navy Blue · Black · Grey
Warranty	7 yr frame · 3 yr electrics · 12 mth footswitch
Certifications	CE · FDA · ISO 13485:2016



 0410 148 292



www.physiomed.au



info@physiomed.au